

# SUPOLENE



SUPOLENE ist ein nicht-resorbierbares, geflochtenes Nahtmaterial aus Polyester, welches zusätzlich mit Teflon (PTFE) beschichtet ist. Durch die hervorragenden Knoteneigenschaften und die dauerhafte Reißfestigkeit kann Supolene überall eingesetzt werden, wo eine dauerhafte Fixierung notwendig ist. Das Nahtmaterial wird in grün und weiß mit einer großen Auswahl von Nadeltypen angeboten.

- nicht-resorbierbarer Polyesterfaden
- extrem zugfest
- hydrophob
- alterungsbeständig
- sehr gewebefreundlich



# SUPOLENE















*SUPOLENE is a non-absorbable multifilament suture based on polyester with additional Teflon (PTFE) coating. Due to the superior knotting properties and the permanent tensile strength Supolene can be used in every procedure where permanent fixation is required. The suture is available in green or white with multiple choices of needle types.*

- *non-absorbable suture*
- *high and permanent tensile strength*
- *water-repellent*
- *provides long-term support in tissue*
- *completely non-irritating*

# SUPOLENE

## Sterile Nadel-Faden-Kombination




Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
<b>1/2-Kreis Rundkörper = HR, HRX</b> <i>1/2 circle, round bodied = HR, HRX</i>						
● HR 12 	5-0	1	0,75	grün/green	9021	36
	4-0	1.5	0,75	grün/green	9022	36
● HR 12 A  Abziehnadel/control release needle	4-0	1.5	8 x 0,45	grün/green	A 90232 Multipack	12
● 2 x HR 15 	5-0	1	0,75	grün/green	9024	36
● HR 17 	5-0	1	0,75	grün/green	9025	36
	4-0	1.5	0,75	grün/green	9026	36
	3-0	2	0,75	grün/green	9027	36
	2-0	3	0,75	grün/green	9028	36
● HR 17 A  Abziehnadel/control release needle	4-0	1.5	8 x 0,45	grün/green	A 90238 Multipack	12
● 2 x HR 17 	4-0	1.5	0,90	grün/green	90210	36
	3-0	2	0,90	grün/green	90214	36
● HR 22 	4-0	1.5	0,75	grün/green	90215	36
	3-0	2	0,75	grün/green	90216	36
	2-0	3	0,75	grün/green	90217	36
● HR 22 A  Abziehnadel/control release needle	4-0	1.5	4 x 0,45	grün/green	A 90233 Multipack	12
● HR 27 	3-0	2	0,75	grün/green	90218	36
	2-0	3	0,75	grün/green	90219	36
	0	3.5	0,75	grün/green	90220	36
	1	4	0,75	grün/green	90221	36
	2-0	3	6 x 0,75	grün/green	M 90222 Multipack	12
	2	5	0,75	grün/green	9041	24
● 2 x HRX 27 	3	6	0,75	grün/green	9042	24
● HR 30 	2-0	3	0,75	grün/green	90227	36
	0	3.5	0,75	grün/green	90228	36
● HR 35 	0	3.5	0,75	grün/green	90230	36
	1	4	0,75	grün/green	90231	36

# SUPOLENE

## Sterile Nadel-Faden-Kombination






Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
<b>1/2-Kreis schneidende Rundkörper = HRT, HRTX</b> <i>1/2 circle, round bodied, trocar point = HRT, HRTX</i>						
<b>▼ HRT 18</b> 	5-0	1	0,45	grün/green	90630	36
	4-0	1.5	0,45	grün/green	90612	36
	5-0	1	0,75	grün/green	9061	36
	4-0	1.5	0,75	grün/green	9062	36
	3-0	2	0,75	grün/green	9063	36
	2-0	3	0,75	grün/green	9064	36
	0	3.5	0,75	grün/green	90652	36
<b>▼ 2 x HRT 18</b> 	4-0	1.5	0,75	grün/green	9067	36
	3-0	2	0,75	grün/green	9068	36
	2-0	3	0,75	grün/green	9069	36
	0	3.5	0,75	grün/green	90610	36
	2-0	3	4 x 0,75	grün/green	M 90614 Multipack	24
	2-0	3	4 x 0,75	grün/weiß green/white	M 90636 Multipack	12
	2-0	3	4 x 0,75	grün/weiß green/white	M 90653 Multipack	12
<b>▼ 2 x HRT 18</b>  + Teflon-Pledgets 6 x 3 mm	2-0	3	8 x 0,75	grün/weiß green/white	M 90656 Multipack	6
	2-0	3	10 x 0,75	grün/weiß green/white	M 90657 Multipack	6
	2-0	3	4 x 0,90	grün/weiß green/white	M 90658 Multipack	12
	0	3.5	4 x 0,75	grün/weiß green/white	M 90659 Multipack	12

# SUPOLENE

## Sterile Nadel-Faden-Kombination










Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces	
	USP	metric					
<b>▼ HRT 25</b> 	3-0	2	0,75	grün/green	90616	36	
	2-0	3	0,75	grün/green	90617	36	
	0	3.5	0,75	grün/green	90618	36	
<b>▼ 2 x HRT 25</b>  + Teflon-Pledgets 6 x 3 mm	2-0	3	4 x 0,75	grün/weiß green/white	M 90654 Multipack	12	
	2-0	3	8 x 0,75	grün/weiß green/white	M 90660 Multipack	6	
	2-0	3	10 x 0,75	grün/weiß green/white	M 90661 Multipack	6	
	2-0	3	4 x 0,90	grün/weiß green/white	M 90662 Multipack	12	
	0	3.5	4 x 0,90	grün/weiß green/white	M 90663 Multipack	12	
	0	3.5	8 x 0,90	grün/weiß green/white	M 90664 Multipack	6	
	<b>▼ HRT 35</b> 	2	5	0,75	grün/green	90655	24
		<b>▼ HRT 45</b> 	2	5	0,75	grün/green	90626
<b>▼ HRTX 45</b>  Cervix-Naht nach MC DONALD Cervix Suture acc. to MC DONALD	6		8	0,75	grün/green	C 90627	12

# SUPOLENE

## Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
<b>▼ HRT 50</b> 	6	8	6 x 0,75	grün/green	M 90631 Multipack	12
<b>▼ HRTX 50</b>  <p>Cervix-Naht nach MC DONALD Cervix Suture acc. to MC DONALD</p>	6	8	0,75	grün/green	C 90628	12
<b>Asymptotische schneidende Rundkörper = ART</b> <i>Asymptotic, round bodied, trocar point = ART</i>						
<b>▼ ART 13</b> 	5-0 4-0	1 1.5	0,45 0,45	grün/green grün/green	9284 9281	12 12
<b>▼ ART 25</b> 	4-0 3-0	1.5 2	0,45 0,45	grün/green grün/green	9282 9283	12 12
<b>3/8-Kreis Rundkörper = DR F</b> <i>3/8 circle, round bodied = DR F</i>						
<b>● DR 10 F</b> 	6-0	0.7	0,75	grün/green	9071	36
<b>● DR 20</b> 	3-0	2	0,75	grün/green	90718	36
<b>Gerade Rundkörper = GR</b> <i>Straight round bodied = GR</i>						
<b>● GR 20</b> 	4-0 3-0	1.5 2	0,45 0,45	grün/green grün/green	9096 9091	36 36
<b>● 2 x GR 20</b> 	4-0	1.5	0,35	grün/green	9093	36
<b>● 2 x GR 65</b> 	3	6	0,75	grün/green	9094	24

# SUPOLENE

## Sterile Nadel-Faden-Kombination











Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
<b>3/8-Kreis schneidende Rundkörper = DRT</b> 3/8 circle, round bodied, trocar point = DRT						
▼ DRT 12 	6-0	0.7	0,75	grün/green	9081	36
	5-0	1	0,75	grün/green	9082	36
▼ DRT 18 	3-0	2	0,75	grün/green	9087	36
	2-0	3	0,75	grün/green	9088	36
<b>1/2-Kreis schneidende = HS</b> 1/2 circle, reverse cutting = HS						
▼ HS 18 	4-0	1.5	0,45	grün/green	9118	36
	3-0	2	0,45	grün/green	9111	36
▼ HS 22 	3-0	2	0,75	grün/green	9112	36
	2-0	3	0,75	grün/green	9113	36
▼ HS 27 	0	3.5	0,75	grün/green	9115	36
▼ HS 35 	1	4	0,75	grün/green	9116	36
▼ HS 40 	1	4	0,75	grün/green	9117	36
<b>3/8-Kreis schneidende = DS, DS F, mit PREMIUM-Schliff = DSM</b> 3/8 circle, reverse cutting = DS, DSF, Premium cut = DSM						
▽ DSM 11  PEREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	grün/green	9151	36

# SUPOLENE

## Sterile Nadel-Faden-Kombination







Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ DS 13 	4-0	1.5	0,45	grün/green	9142	36
▽ DSM 13  PEREMIUM-Schliff / PREMIUM cut	5-0 4-0	1 1.5	0,45 0,45	grün/green grün/green	9152 9153	36 36
▽ DSM 16  PEREMIUM-Schliff / PREMIUM cut	5-0 4-0	1 1.5	0,45 0,45	grün/green grün/green	9155 9156	36 36
▼ DS 18 	5-0 4-0 3-0	1 1.5 2	0,45 0,45 0,45	grün/green grün/green grün/green	9144 9145 9146	36 36 36
▼ DS 18 F 	4-0	1.5	0,45	grün/green	91427	36
▽ DSM 18 	4-0 3-0 2-0 4-0 3-0	1.5 2 3 1.5 2	0,45 0,45 0,45 0,75 4 x 0,45	grün/green grün/green grün/green grün/green grün/green	9157 9158 9159 91526 M 91515 Multipack	36 36 36 36 24
▼ DS 21 	4-0 3-0 2-0	1.5 2 3	0,45 0,45 0,45	grün/green grün/green grün/green	91423 91424 91425	36 36 36
▽ DSM 21 	4-0 3-0 2-0	1.5 2 3	0,45 0,45 0,45	grün/green grün/green grün/green	91516 91517 91518	36 36 36
▼ DS 24 	4-0 3-0 2-0 3-0 0	1.5 2 3 2 3.5	0,45 0,45 0,45 0,75 0,75	grün/green grün/green grün/green grün/green grün/green	9147 9148 9149 91410 91433	36 36 36 36 36
▽ DSM 24  PEREMIUM-Schliff / PREMIUM cut	3-0 2-0	2 3	0,75 0,75	grün/green grün/green	91513 91514	36 36


# SUPOLENE

## Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ DS 27 	2-0	3	0,45	grün/green	91411	36
	2-0	3	0,75	grün/green	91426	36
▽ DSM 27  PEREMIUM-Schliff / PREMIUM cut	2-0	3	0,45	grün/green	91520	36
	2-0	3	0,75	grün/green	91519	36
▼ DS 30 	2-0	3	0,75	grün/green	91412	36
	0	3.5	0,75	grün/green	91413	36
	1	4	0,75	grün/green	91428	36
▽ DSM 30  PEREMIUM-Schliff / PREMIUM cut	2-0	3	0,75	grün/green	91521	36
	0	3.5	0,75	grün/green	91522	36
▼ DS 40 	2-0	3	0,75	grün/green	91435	36
	0	3.5	0,75	grün/green	91414	36
▼ DS 50 	1	4	0,75	grün/green	91434	36

### ▼ 2 x DS 90



Doppelte Unterstützungsnaht  
Tension Suture, double-armed

Zubehör: Je 2 Nylonplatten 80 x 23 mm  
und 2 Moosgummiplatten 82 x 25 mm  
Incl. 4 buttons

5 7  
6 8

0,75  
0,75

grün/green  
grün/green

91415 U  
91416 U





6  
6



# SUPOLENE

## Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
<b>Gerade schneidende = GS</b> <i>Straight, cutting = GS</i>						
▼ <b>GS 60</b> 	0	3.5	0,45	grün/green	9178	36
	2-0	3	0,75	grün/green	9174	36
<b>1/2-Kreis gebogene schneidende = KS</b> <i>1/2 circle, curved, reverse cutting = KS</i>						
▼ <b>KS 20</b> 	3-0	2	0,45	grün/green	9181	36
<b>1/2-Kreis schneidende Spatula, Premium-Schliff = HSPM, HSPM F</b> <i>1/2 circle, spatulated, Premium Cut = HSPM, HSPM F</i>						
▽ <b>2 x HSPM 8</b>  PEREMIUM-Schliff / <i>PREMIUM cut</i>	4-0	1.5	0,45	weiß/white	9191	24
<b>1/4-Kreis schneidende Spatula, Premium-Schliff = VSPM F</b> <i>1/4 circle, spatulated, Premium Cut = VSPM F</i>						
▽ <b>2 x VSPM 8 F</b>  PEREMIUM-Schliff / <i>PREMIUM cut</i>	6-0	0.7	0,45	weiß/white	9211	24
	5-0	1	0,45	weiß/white	9212	24

# SUPOLENE

Einzelfadenpackungen, steril in vorgeschnittenen Längen

Needle - suture - combination, sterile

Stärke Nr. gauge no. USP    metric		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
3-0	2	5 x 0,35	grün/green	W 453	24
2-0	3	5 x 0,35	grün/green	W 553	24
0	3.5	5 x 0,35	grün/green	W 653	24
4-0	1.5	3 x 0,50	grün/green	W 335	24
3-0	2	3 x 0,50	grün/green	W 435	24
2-0	3	3 x 0,50	grün/green	W 535	24
0	3.5	3 x 0,50	grün/green	W 635	24
1	4	3 x 0,50	grün/green	W 735	24
2	5	3 x 0,50	grün/green	W 835	24
4-0	1.5	5 x 0,50	grün/green	W 355	24
3-0	2	5 x 0,50	grün/green	W 455	24
2-0	3	5 x 0,50	grün/green	W 555	24
0	3.5	5 x 0,50	grün/green	W 655	24
1	4	5 x 0,50	grün/green	W 755	24
2	5	5 x 0,50	grün/green	W 855	24
4-0	1.5	6 x 0,50	grün/green	SW 355	24
3-0	2	6 x 0,50	grün/green	SW 455	24
2-0	3	6 x 0,50	grün/green	SW 555	24
0	3.5	6 x 0,50	grün/green	SW 655	24
1	4	6 x 0,50	grün/green	SW 755	24
2	5	6 x 0,50	grün/green	SW 855	24
4-0	1.5	10 x 0,50	grün/green	W 3105	24
3-0	2	10 x 0,50	grün/green	W 4105	24
2-0	3	10 x 0,50	grün/green	W 5105	24
2-0	3	12 x 0,50	grün/green	SW 5105	24
0	3.5	12 x 0,50	grün/green	SW 6105	24
1	4	12 x 0,50	grün/green	SW 7105	24
3-0	2	17 x 0,50	grün/green	W 4175	24
2-0	3	2 x 0,75	grün/green	W 520	24
0	3.5	2 x 0,75	grün/green	W 620	24
1	4	2 x 0,75	grün/green	W 720	24
2	5	2 x 0,75	grün/green	W 820	24
3	6	2 x 0,75	grün/green	W 920	24
1	4	5 x 0,75	grün/green	W 750	24
2	5	5 x 0,75	grün/green	W 850	24
3	6	5 x 0,75	grün/green	W 950	24