

SUPRAMID



SUPRAMID ist ein nicht-resorbierbares Nahtmaterial aus Polyamid. Ursprünglich ein monofiler Faden, werden die Faserbündel ab der USP 4-0 und dicker von einem Polyamid-Mantel umhüllt. Dadurch erreicht das „Pseudomonofilament“ Handlungseigenschaften wie ein geflochtenes Nahtmaterial mit dem geringeren Gewebetrauma eines Monofilaments. Das Nahtmaterial wird in schwarz und weiß angeboten.

- nicht-resorbierbarer Polyamid-Faden
- sehr geschmeidig
- zugfest
- geringe Gewebereaktion
- gute Knoteneigenschaften



SUPRAMID











SUPRAMID is a non-absorbable suture based on polyamid. Basically a monofilament the single suture threads in USP 4-0 and thicker are wrapped with polyamid. Through this the "pseudomonofilament" can be used as a multifilament with the minimized tissue drag comparable to a monofilament.

- *non-absorbable polyamide-suture*
- *very soft*
- *high tensile strength*
- *less tissue reaction*
- *good knotting properties*

SUPRAMID

Sterile Nadel-Faden-Kombination





Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
1/2-Kreis schneidende Rundkörper = HRT <i>1/2 circle, round bodied, trocar point = HRT</i>						
▼ HRT 18 	4-0	1.5	0,50	schwarz/black	6061	36
	3-0	2	0,50	schwarz/black	6062	36
	4-0	1.5	0,75	schwarz/black	6063	36
	3-0	2	0,75	schwarz/black	6064	36
Gerade Rundkörper = GR <i>straight, round bodied = GR</i>						
● GR 20 	5-0	1.5	0,45	schwarz/black	6091	36
	4-0	1	0,45	schwarz/black	6090	36
1/2-Kreis schneidende = HS <i>1/2 circle, reverse cutting = HS</i>						
▼ HS 15 	5-0	1	0,45	schwarz/black	6119	36
	4-0	1.5	0,45	schwarz/black	6118	36
▼ HS 18 	5-0	1	0,45	schwarz/black	6111	36
	4-0	1.5	0,45	schwarz/black	6112	36
	3-0	2	0,45	schwarz/black	6114	36
	4-0	1.5	0,75	schwarz/black	61111	36
	3-0	2	0,75	schwarz/black	6117	36
▼ HS 20 	4-0	1.5	0,75	schwarz/black	61110	36
	3-0	2	0,75	schwarz/black	6116	36
▼ HS 22 	3-0	2	0,45	schwarz/black	6113	36
	2-0	3	0,75	schwarz/black	6115	36
3/8-Kreis Rundkörper = DR <i>3/8 circle, round bodied = DR</i>						
● DR 18 / Schlinge / loop 	4-0	1.5	0,50	weiß/white	6070	36
Tsuge-Naht / Tsuge Suture						
3/8-Kreis schneidende Rundkörper = DRT <i>3/8 circle, trocar point, round bodied = DRT</i>						
▼ DRT 18 / Schlinge / loop 	4-0	1.5	0,50	weiß/white	6080	36
Tsuge-Naht / Tsuge Suture						

SUPRAMID

Sterile Nadel-Faden-Kombination











Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
3/8-Kreis schneidende = DS, DS F mit PREMIUM-Schliff = DSM <i>3/8 circle, reverse cutting = DS, DS F PREMIUM cut = DSM</i>						
 DSM 11 PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz/black	6151	36
	5-0	1	0,45	schwarz/black	6152	36
 DS 13	5-0	1	0,45	schwarz/black	6141	36
 DSM 13 PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	schwarz/black	6158	36
	4-0	1.5	0,45	schwarz/black	6154	36
 DSM 16 PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	schwarz/black	6157	36
	4-0	1.5	0,75	schwarz/black	6153	36

SUPRAMID

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
	USP	metric				
▼ DS 18 	4-0	1.5	0,45	schwarz/black	6142	36
	3-0	2	0,45	schwarz/black	6143	36
	4-0	1.5	0,75	schwarz/black	61428	36
	3-0	2	0,75	schwarz/black	61429	36
	3-0	2	6 x 0,45	schwarz/black	M 6144 Multipack	12
	4-0	1.5	0,45	weiß/white	61430	36
	3-0	2	0,45	weiß/white	61431	36
▼ DS 18 F 	4-0	1.5	0,45	schwarz/black	61420	36
▽ DSM 18 	4-0	1.5	0,75	schwarz/black	6155	36
	3	2	0,75	schwarz/black	6156	36
▼ DS 24 	4-0	1.5	0,45	schwarz/black	6146	36
	3-0	2	0,45	schwarz/black	6147	36
	2-0	3	0,45	schwarz/black	61422	36
	3-0	2	0,75	schwarz/black	61425	36
	3-0	2	6 x 0,45	schwarz/black	M 6148 Multipack	12
▼ DS 24 F 	5-0	1	0,45	schwarz/black	61418	36
▼ DS 30 	2-0	3	0,45	schwarz/black	6149	36
	0	3.5	0,45	schwarz/black	61423	36
	2-0	3	0,75	schwarz/black	61426	36
	0	3.5	0,75	schwarz/black	61427	36
	2-0	3	6 x 0,45	schwarz/black	M 61411 Multipack	12
▼ DS 40 	3-0	2	0,45	schwarz/black	61412	36
	2-0	2.5	0,45	schwarz/black	61413	36
	0	3.5	0,75	weiß/white	61421	36
<p>1/4-Kreis schneidende Spatula, PREMIUM-Schliff = VSPM, VSPM F <i>1/4 circle, spatulated, PREMIUM cut = VSPM, VSPM F</i></p>						
▽ VSPM 8  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,75	weiß/white	6211	24
▽ 2x VSPM 8  PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,25	weiß/white	6213	24
	4-0	1.5	0,45	weiß/white	6214	24
▽ 2 x VSPM 8 F  PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	weiß/white	6212	24

SUPRAMID

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-neededled

Stärke Nr. gauge no. USP metric		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Stück/pieces
4-0	1.5	3 x 0,50	schwarz/ <i>black</i>	SAM 335	24
3-0	2	3 x 0,50	schwarz/ <i>black</i>	SAM 435	24
2-0	3	3 x 0,50	schwarz/ <i>black</i>	SAM 535	24
0	3.5	3 x 0,50	schwarz/ <i>black</i>	SAM 635	24
1	4	3 x 0,50	schwarz/ <i>black</i>	SAM 735	24
4-0	1.5	3 x 0,50	weiß/ <i>white</i>	M 335	24
3-0	2	3 x 0,50	weiß/ <i>white</i>	M 435	24
2-0	3	3 x 0,50	weiß/ <i>white</i>	M 535	24
0	3.5	3 x 0,50	weiß/ <i>white</i>	M 635	24
1	4	3 x 0,50	weiß/ <i>white</i>	M 735	24
4-0	1.5	5 x 0,50	schwarz/ <i>black</i>	SAM 355	24
3-0	2	5 x 0,50	schwarz/ <i>black</i>	SAM 455	24
2-0	3	5 x 0,50	schwarz/ <i>black</i>	SAM 555	24
0	3.5	5 x 0,50	schwarz/ <i>black</i>	SAM 655	24
3-0	2	5 x 0,50	weiß/ <i>white</i>	M 455	24
2-0	3	5 x 0,50	weiß/ <i>white</i>	M 555	24
0	3.5	5 x 0,50	weiß/ <i>white</i>	M 655	24
3-0	2	10 x 0,50	schwarz/ <i>black</i>	SAM 4105	24
3-0	2	10 x 0,50	weiß/ <i>white</i>	M 4105	24
2-0	3	10 x 0,50	weiß/ <i>white</i>	M 5105	24